




IVTScrip™ mRNA-Human APOE, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25231MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25231MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the APOE protein. This product can be used in Astrocyte-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000041.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000041.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	AD2; LPG; APO-E; ApoE4; LDLCQ5
Description	The protein encoded by this gene is a major apoprotein of the chylomicron. It binds to a specific liver and peripheral cell receptor, and is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. This gene maps to chromosome 19 in a cluster with the related apolipoprotein C1 and C2 genes. Mutations in this gene result in familial dysbetalipoproteinemia, or type III hyperlipoproteinemia (HLP III), in which increased plasma cholesterol and triglycerides are the consequence of impaired clearance of chylomicron and VLDL remnants. [provided by RefSeq, Jun 2016]