




IVTScrip™ mRNA-Human ARHGEF38, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25318MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25318MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ARHGEF38 protein. This product can be used in Ciliated epithelial cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001242729.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001242729.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARHGEF38 (Rho Guanine Nucleotide Exchange Factor 38) is a Protein Coding gene. Among its related pathways are Signaling by GPCR and p75 NTR receptor-mediated signalling. Gene Ontology (GO) annotations related to this gene include Rho guanyl-nucleotide exchange factor activity. An important paralog of this gene is DNMBP.
--------------------	---