



IVTScrip™ mRNA-Human ARRDC4, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25382MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25382MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ARRDC4 protein. This product can be used in Alpha cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_183376.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_183376.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARRDC4 (Arrestin Domain Containing 4) is a Protein Coding gene. Diseases associated with ARRDC4 include Mitochondrial Myopathy, Encephalopathy, Lactic Acidosis, And Stroke-Like Episodes. An important paralog of this gene is ARRDC3.
--------------------	---