



IVTScrip™ mRNA-Human ASRGL1, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25427MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25427MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ASRGL1 protein. This product can be used in Mitotic fetal germ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001083926.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001083926.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ALP; ALP1; CRASH
Description	ASRGL1 (Asparaginase And Isoaspartyl Peptidase 1) is a Protein Coding gene. Diseases associated with ASRGL1 include Aspartylglucosaminuria and Immunodeficiency 48. Among its related pathways are Histidine, lysine, phenylalanine, tyrosine, proline and tryptophan catabolism and Viral mRNA Translation. Gene Ontology (GO) annotations related to this gene include hydrolase activity and asparaginase activity. An important paralog of this gene is TASP1.