



IVTScrip™ mRNA-Human ATP1A1-AS1, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25483MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25483MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ATP1A1-AS1 protein. This product can be used in FGFR1HighNME5- epithelial cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ATP1A1OS; C1orf203
Description	ATP1A1-AS1 (ATP1A1 Antisense RNA 1) is an RNA Gene, and is affiliated with the lncRNA class. Diseases associated with ATP1A1-AS1 include Charcot-Marie-Tooth Disease, Axonal, Type 2Dd and Hypomagnesemia, Seizures, And Mental Retardation 2.