



IVTScrip™ mRNA-Human ARL6IP4, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25648MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25648MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the ARL6IP4 protein. This product can be used in Retinoid acid signaling-responsive fetal germ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001002251.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001002251.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	SR-25; SRp25; SFRS20; SRp37
Description	ARL6IP4 (ADP Ribosylation Factor Like GTPase 6 Interacting Protein 4) is a Protein Coding gene. Diseases associated with ARL6IP4 include Tabes Dorsalis and Tricuspid Valve Stenosis.