



IVTScrip™ mRNA-Human ARPC5L, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25672MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25672MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the ARPC5L protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_030978.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_030978.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ARC16-2
Description	ARPC5L (Actin Related Protein 2/3 Complex Subunit 5 Like) is a Protein Coding gene. Diseases associated with ARPC5L include Kleptomania and Leopard Syndrome 1. Among its related pathways are Yersinia infection and Shigellosis. Gene Ontology (GO) annotations related to this gene include actin binding. An important paralog of this gene is ARPC5.