



IVTScrip™ mRNA-Human ASXL2, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25732MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25732MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the ASXL2 protein. This product can be used in Regulatory T (Treg) cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001369346.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001369346.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ASXH2; SHAPNS
Description	This gene encodes a member of a family of epigenetic regulators that bind various histone-modifying enzymes and are involved in the assembly of transcription factors at specific genomic loci. Naturally occurring mutations in this gene are associated with cancer in several tissue types (breast, bladder, pancreas, ovary, prostate, and blood). This gene plays an important role in neurodevelopment, cardiac function, adipogenesis, and osteoclastogenesis. [provided by RefSeq, Feb 2017]