



IVTScrip™ mRNA-Human ARHGAP20, (Cap 1, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK25886MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25886MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the ARHGAP20 protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001258415.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001258415.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	RARHOGAP
Description	The protein encoded by this gene is an activator of RHO-type GTPases, transducing a signal from RAP1 to RHO and impacting neurite outgrowth. [provided by RefSeq, Sep 2016]