




IVTScrip™ mRNA-Human CD56, (Cap 0, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2610MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2610MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the CD56 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000615.6
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000615.6
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	NCAM; MSK39; NCAM-1; NCAM; CD56
Description	This gene is a protein coding gene. Exhibits identical protein binding activity. Predicted to be involved in several processes, including commissural neuron axon guidance; regulation of exocyst assembly; and regulation of signal transduction. Localizes to cytosol and external side of plasma membrane. Implicated in bipolar disorder and pancreatic cancer. Biomarker of adenoma; cholangiocarcinoma; insulinoma; pancreatic cancer; and type 1 diabetes mellitus.