




IVTScrip™ mRNA-Human ARHGAP10, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK26181MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK26181MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the ARHGAP10 protein. This product can be used in FGFR1 High NME5- epithelial cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024605.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024605.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARHGAP10 (Rho GTPase Activating Protein 10) is a Protein Coding gene. Diseases associated with ARHGAP10 include Epidermolysis Bullosa Simplex With Nail Dystrophy. Among its related pathways are Signaling by GPCR and p75 NTR receptor-mediated signalling. Gene Ontology (GO) annotations related to this gene include GTPase activator activity. An important paralog of this gene is ARHGAP26.
--------------------	---