




IVTScrip™ mRNA-Human ATCAY, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK26343MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK26343MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the ATCAY protein. This product can be used in Oligodendrocyte-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_033064.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_033064.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	CLAC; BNIP-H
Description	This gene encodes a neuron-restricted protein that contains a CRAL-TRIO motif common to proteins that bind small lipophilic molecules. Mutations in this gene are associated with cerebellar ataxia, Cayman type. [provided by RefSeq, Jul 2008]