




IVTScrip™ mRNA-Human ARHGAP45, (Cap 1, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK26500MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK26500MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the ARHGAP45 protein. This product can be used in Naive CD8+ T cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001258328.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001258328.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARHGAP45 (Rho GTPase Activating Protein 45) is a Protein Coding gene. Diseases associated with ARHGAP45 include Graft-Versus-Host Disease and Acute Graft Versus Host Disease. Among its related pathways are Signaling by GPCR and p75 NTR receptor-mediated signalling. An important paralog of this gene is ARHGAP29.
--------------------	--