



IVTScrip™ mRNA-Human ARHGAP30, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK27095MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK27095MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ARHGAP30 protein. This product can be used in FGFR1 High NME5- epithelial cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001025598.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001025598.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARHGAP30 (Rho GTPase Activating Protein 30) is a Protein Coding gene. Among its related pathways are Signaling by GPCR and p75 NTR receptor-mediated signalling. Gene Ontology (GO) annotations related to this gene include GTPase activator activity. An important paralog of this gene is ARHGAP32.
--------------------	--