



## IVTScrip™ mRNA-Human ATF2, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK27245MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK27245MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ATF2 protein. This product can be used in Natural killer T (NKT) cell-related researches.

#### Specifications

<b>Modified bases</b>	Pseudo-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001256090.2
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

#### SPECIFICATIONS

<b>Modified bases</b>	Pseudo-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001256090.2
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

## GENE INFORMATION

<b>Alternative Names</b>	HB16; CREB2; TREB7; CREB-2; CRE-BP1
<b>Description</b>	<p>This gene encodes a transcription factor that is a member of the leucine zipper family of DNA binding proteins. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. This protein binds to the cAMP-responsive element (CRE), an octameric palindrome. It forms a homodimer or a heterodimer with c-Jun and stimulates CRE-dependent transcription. This protein is also a histone acetyltransferase (HAT) that specifically acetylates histones H2B and H4 in vitro; thus it may represent a class of sequence-specific factors that activate transcription by direct effects on chromatin components. The encoded protein may also be involved in cell's DNA damage response independent of its role in transcriptional regulation. Several alternatively spliced transcript variants have been found for this gene [provided by RefSeq, Jan 2014]</p>