



IVTScrip™ mRNA-Human ATF7IP, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK27249MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK27249MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ATF7IP protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001286514.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001286514.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	AM; MCAF; p621; MCAF1; ATF-IP
Description	ATF7IP is a multifunctional nuclear protein that associates with heterochromatin. It can act as a transcriptional coactivator or corepressor depending upon its binding partners (summary by Liu et al., 2009 [PubMed 19106100]).[supplied by OMIM, Nov 2010]