




IVTScrip™ mRNA-Human ATG101, (Cap 1, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK27852MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK27852MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the ATG101 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001098673.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001098673.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	C12orf44
Description	ATG101 (Autophagy Related 101) is a Protein Coding gene. Diseases associated with ATG101 include Huntington Disease. Among its related pathways are Pathways of neurodegeneration - multiple diseases and Macroautophagy. Gene Ontology (GO) annotations related to this gene include identical protein binding.