



IVTScrip™ mRNA-Human ATAD5, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK28141MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK28141MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the ATAD5 protein. This product can be used in Goblet cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024857.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024857.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ELG1; FRAG1; C17orf41
Description	ATAD5 (ATPase Family AAA Domain Containing 5) is a Protein Coding gene. Diseases associated with ATAD5 include Chromosome 17Q11.2 Deletion Syndrome, 1.4-Mb and Cafe-Au-Lait Spots, Multiple. An important paralog of this gene is CRLF3.