



## IVTScrip™ mRNA-Human ARF4-AS1, (Cap 1, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK28266MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK28266MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the ARF4-AS1 protein. This product can be used in Natural killer T (NKT) cell-related researches.

#### Specifications

<b>Modified bases</b>	5-Methyl-CTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

#### SPECIFICATIONS

<b>Modified bases</b>	5-Methyl-CTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

## GENE INFORMATION

**Description** ARF4-AS1 (ARF4 Antisense RNA 1) is an RNA Gene, and is affiliated with the lncRNA class.