



## IVTScrip™ mRNA-Human APP, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK28541MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK28541MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the APP protein. This product can be used in Beta cell-related researches.

#### Specifications

|                |                           |
|----------------|---------------------------|
| Modified bases | 5-Methyl-CTP & Pseudo-UTP |
| 5' Cap         | Cap 1                     |
| Species        | Human                     |
| RefSeq         | NM_000484.4               |
| Applications   | Gene therapy research     |
| Format         | Powder                    |
| Quantity       | 100 µg                    |
| Purification   | Chromatography            |

#### SPECIFICATIONS

|                       |                           |
|-----------------------|---------------------------|
| <b>Modified bases</b> | 5-Methyl-CTP & Pseudo-UTP |
| <b>5' Cap</b>         | Cap 1                     |
| <b>Species</b>        | Human                     |
| <b>RefSeq</b>         | NM_000484.4               |
| <b>Applications</b>   | Gene therapy research     |
| <b>Format</b>         | Powder                    |
| <b>Quantity</b>       | 100 µg                    |
| <b>Purification</b>   | Chromatography            |

## GENE INFORMATION

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Description</b>       | <p>This gene encodes a cell surface receptor and transmembrane precursor protein that is cleaved by secretases to form a number of peptides. Some of these peptides are secreted and can bind to the acetyltransferase complex APBB1/TIP60 to promote transcriptional activation, while others form the protein basis of the amyloid plaques found in the brains of patients with Alzheimer disease. In addition, two of the peptides are antimicrobial peptides, having been shown to have bacteriocidal and antifungal activities. Mutations in this gene have been implicated in autosomal dominant Alzheimer disease and cerebroarterial amyloidosis (cerebral amyloid angiopathy). Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2014]</p> |