



IVTScrip™ mRNA-Human ASB13, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK28695MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK28695MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the ASB13 protein. This product can be used in B cell-related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024701.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024701.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and a SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Multiple alternatively spliced transcript variants, both protein-coding and not protein-coding, have been described for this gene. [provided by RefSeq, Nov 2010]
--------------------	---