



IVTScrip™ mRNA-Human CCR3, (Cap 0, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2972MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2972MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the CCR3 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_178329
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_178329
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CD193; CMKBR3; CCR-3
Description	This gene is a protein coding gene. Exhibits C-C chemokine binding activity. Involved in positive regulation of angiogenesis. Predicted to localize to cytoplasm and external side of plasma membrane. Implicated in asthma; ovary epithelial cancer; and rhinitis. Biomarker of several diseases, including autoimmune disease (multiple); cystic fibrosis; hepatitis B; lung disease (multiple); and periaapical granuloma.