



## IVTScip™ mRNA-Human ALDH1, (Cap 0, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2999MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK2999MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ALDH1 protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	2-Thio-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_000689
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

### SPECIFICATIONS

<b>Modified bases</b>	2-Thio-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_000689
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

## GENE INFORMATION

<b>Alternative Names</b>	ALDC; ALDH11; ALDH-E1
<b>Description</b>	This gene is a protein coding gene. Exhibits aldehyde dehydrogenase (NAD <sup>+</sup> ) activity. Involved in oxidation-reduction process. Localizes to cytosol. Biomarker of pancreatic cancer.