




IVTScrip™ mRNA-Human D105, (Cap 0, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK3128MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3128MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the D105 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000118
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000118
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	ORW1; CD105; ORW; FLJ41744; HHT1; END
Description	This gene is a protein coding gene. Exhibits several functions, including galactose binding activity; protein homodimerization activity; and transforming growth factor beta receptor binding activity. Contributes to BMP binding activity. Involved in several processes, including central nervous system vasculogenesis; regulation of transcription by RNA polymerase II; and regulation of transmembrane receptor protein serine/threonine kinase signaling pathway. Localizes to several cellular components, including external side of plasma membrane; extracellular space; and nuclear body. Implicated in arteriovenous malformation; arteriovenous malformations of the brain; breast cancer; hereditary hemorrhagic telangiectasia; and intracranial aneurysm. Biomarker of several diseases, including artery disease (multiple); clear cell renal cell carcinoma; diabetic retinopathy; hematologic cancer (multiple); and secondary hyperparathyroidism.