



IVTScrip™ mRNA-Human FOXP3, (Cap 0, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK3188MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3188MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the FOXP3 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_014009
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_014009
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	XPID; DIETER; PIDX
Description	This gene is a protein coding gene. Exhibits several functions, including NFAT protein binding activity; enzyme binding activity; and protein homodimerization activity. Involved in several processes, including negative regulation of T cell activation; negative regulation of cytokine production; and regulation of transcription, DNA-templated. Localizes to cytoplasm and nucleoplasm. Implicated in colorectal cancer; immune dysregulation-polyendocrinopathy-enteropathy-X-linked syndrome; pulmonary tuberculosis; and type 1 diabetes mellitus. Biomarker of several diseases, including adult T-cell leukemia; chlamydia; cutaneous leishmaniasis; liver disease (multiple); and tuberculosis (multiple).