




IVTScrip™ mRNA-Human CD24, (Cap 0, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK3227MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3227MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the CD24 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001291737
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001291737
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MGC75043; CD24A
Description	This gene is a protein coding gene. Exhibits protein kinase binding activity and protein tyrosine kinase activator activity. Involved in several processes, including glomerular epithelial cell differentiation; positive regulation of T cell activation; and positive regulation of nephron tubule epithelial cell differentiation. Localizes to cell surface; intracellular membrane-bounded organelle; and membrane raft. Implicated in multiple sclerosis.