



## IVTScrip™ mRNA-Human CD6, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK3485MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK3485MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CD6 protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	5-Methyl-CTP & Pseudo-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001254750
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

### SPECIFICATIONS

<b>Modified bases</b>	5-Methyl-CTP & Pseudo-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001254750
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

## GENE INFORMATION

<b>Alternative Names</b>	FLJ44171; TP120
<b>Description</b>	This gene is a protein coding gene. Exhibits identical protein binding activity; lipopolysaccharide binding activity; and lipoteichoic acid binding activity. Involved in several processes, including heterophilic cell-cell adhesion via plasma membrane cell adhesion molecules; lipopolysaccharide-mediated signaling pathway; and positive regulation of cytokine production involved in inflammatory response. Localizes to immunological synapse and intrinsic component of plasma membrane. Colocalizes with T cell receptor complex.