




IVTScrip™ mRNA-Human LIN2, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK3563MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3563MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the LIN2 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001126054
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001126054
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	LIN2; MRXSNA; FLJ22219; MICPCH
Description	This gene is a protein coding gene. Exhibits protein serine/threonine kinase activity. Involved in several processes, including negative regulation of cell-matrix adhesion; negative regulation of keratinocyte proliferation; and negative regulation of wound healing. Localizes to several cellular components, including basement membrane; nuclear periphery; and nucleolus. Implicated in FG syndrome and syndromic X-linked intellectual disability Najm type.