



IVTScrip™ mRNA-Human CD203c, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK3573MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3573MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CD203c protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_005021
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_005021
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CD203c; PDNP3; NPP3
Description	This gene is a protein coding gene. Exhibits metal ion binding activity; nucleoside-triphosphate diphosphatase activity; and phosphodiesterase I activity. Involved in ATP metabolic process; nucleoside triphosphate catabolic process; and pyrimidine nucleotide metabolic process. Localizes to external side of plasma membrane and perinuclear region of cytoplasm.