



IVTScrip™ mRNA-Human IL1RAP, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK3598MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3598MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the IL1RAP protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002182
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002182
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	C3orf13; IL-1R3; IL1R3
Description	This gene is a protein coding gene. Predicted to have interleukin-1 receptor activity; interleukin-1 receptor binding activity; and protein tyrosine kinase binding activity. Predicted to be involved in positive regulation of NF-kappaB transcription factor activity; positive regulation of dendrite development; and positive regulation of synapse assembly. Localizes to membrane.