



IVTScrip™ mRNA-Human MEPE, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK3616MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3616MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the MEPE protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001184694
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001184694
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	OF45
Description	This gene is a protein coding gene. Predicted to have extracellular matrix protein binding activity. Predicted to be involved in biomineral tissue development and negative regulation of bone mineralization. Predicted to localize to extracellular matrix.