



IVTScrip™ mRNA-Human CD49d, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK3876MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK3876MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the CD49d protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000885
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000885
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CD49D; IA4
Description	This gene is a protein coding gene. Exhibits cell adhesion molecule binding activity. Contributes to C-X3-C chemokine binding activity. Involved in several processes, including cell-cell adhesion in response to extracellular stimulus; cell-cell adhesion mediated by integrin; and cellular response to a myloid-beta. Localizes to cell surface and integrin alpha4-beta7 complex. Implicated in autism spectrum disorder and pancreatic cancer.