



IVTScrip™ mRNA-Human CD9, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4035MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4035MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It codes the CD9 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001330312.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001330312.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MIC3
Description	This gene is a protein coding gene. Exhibits integrin binding activity. Involved in several processes, including fusion of sperm to egg plasma membrane involved in single fertilization; glial cell migration; and negative regulation of platelet aggregation. Localizes to extracellular exosome and protein-containing complex. Implicated in prostate adenocarcinoma. Biomarker of breast carcinoma; female reproductive organ cancer (multiple); severe acute respiratory syndrome; transitional cell carcinoma; and urinary bladder cancer.