




IVTScrip™ mRNA-Human CDH11, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4135MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4135MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It codes the CDH11 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001308392
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001308392
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	OB; OSF-4; ESWS
Description	This gene is a protein coding gene. Predicted to have cadherin binding activity and calcium ion binding activity. Predicted to be involved in several processes, including calcium-dependent cell-cell adhesion via plasma membrane cell adhesion molecules; cell morphogenesis; and cell-cell junction organization. Localizes to extracellular exosome. Implicated in Elsahy-Waters syndrome.