



IVTScrip™ mRNA-Human CD62L, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4157MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4157MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It codes the CD62L protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000655
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000655
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	LECAM1; PLNHR; LEU8; LAM1
Description	This gene is a protein coding gene. Exhibits calcium ion binding activity; oligosaccharide binding activity; and protease binding activity. Involved in calcium-dependent cell-cell adhesion via plasma membrane cell adhesion molecules and leukocyte tethering or rolling. Localizes to integral component of plasma membrane. Implicated in Crohn's disease; IgA glomerulonephritis; and ulcerative colitis. Biomarker of Alzheimer's disease; IgA glomerulonephritis; and systemic lupus erythematosus.