



IVTScrip™ mRNA-Human CCR6, (Cap 1, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4329MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4329MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the CCR6 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_031409
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_031409
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CMKBR6; DRY6
Description	This gene is a protein coding gene. Exhibits C-C chemokine binding activity and C-C chemokine receptor activity. Involved in calcium-mediated signaling; cell chemotaxis; and positive regulation of flagellated sperm motility involved in capacitation. Localizes to cell surface; sperm flagellum; and sperm plasma membrane. Biomarker of several diseases, including IgA glomerulonephritis; cutaneous T cell lymphoma; exanthem; psoriasis; and systemic scleroderma.