



IVTScrip™ mRNA-Human PDGFC, (Cap 1, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4358MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4358MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the PDGFC protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_016205
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_016205
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	SCDGF; VEGF-E
Description	This gene is a protein coding gene. Exhibits platelet-derived growth factor receptor binding activity and protein homodimerization activity. Involved in platelet-derived growth factor receptor signaling pathway and positive regulation of fibroblast proliferation. Localizes to several cellular components, including cell surface; cytosol; and extracellular space.