



IVTScrip™ mRNA-Human AID, (Cap 1, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4444MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4444MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the AID protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_020661
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_020661
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	ARP2
Description	This gene is a protein coding gene. Exhibits cytidine deaminase activity; identical protein binding activity; and ubiquitin protein ligase binding activity. Involved in DNA demethylation; regulation of cellular macromolecule biosynthetic process; and somatic hypermutation of immunoglobulin genes. Localizes to cytoplasm; nucleus; and protein-containing complex. Implicated in hyperimmunoglobulin syndrome and immunodeficiency with hyper-IgM type 2. Biomarker of acute lymphoblastic leukemia and swine influenza.