




IVTScrip™ mRNA-Human CD127, (Cap 1, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4450MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4450MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the CD127 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_002185
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_002185
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	IL-7RA; CD127; CDW127; IL7RA
Description	This gene is a protein coding gene. Exhibits cytokine receptor activity. Involved in positive regulation of cell population proliferation and positive regulation of receptor signaling pathway via STAT. Predicted to localize to clathrin-coated vesicle membrane and plasma membrane. Implicated in severe combined immunodeficiency and severe combined immunodeficiency, autosomal recessive, T cell-negative, B cell-positive, Nk cell-positive.