



IVTScrip™ mRNA-Human CD147, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4595MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4595MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the CD147 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_0011741057
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_0011741057
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	SLC7A11; EMPRIN
Description	This gene is a protein coding gene. Exhibits signaling receptor activity and virus receptor activity. Involved in several processes, including neutrophil chemotaxis; photoreceptor cell maintenance; and positive regulation of cytokine production. Localizes to several cellular components, including endoplasmic reticulum membrane; endosome; and plasma membrane. Biomarker of breast ductal carcinoma; carcinoma (multiple); and cervical cancer.