




IVTScrip™ mRNA-Human CD1d, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4604MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4604MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the CD1d protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001766
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001766
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	R3G1; R3
Description	This gene is a protein coding gene. Exhibits cell adhesion molecule binding activity; exogenous lipid antigen binding activity; and histone binding activity. Involved in antigen processing and presentation, endogenous lipid antigen via MHC class Ib and positive regulation of T cell proliferation. Localizes to cell surface and endoplasmic reticulum.