




IVTScrip™ mRNA-Human CD71, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4748MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4748MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the CD71 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001313965
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001313965
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	T9; TR; CD71; TFR; TFR1; TRFR
Description	This gene is a protein coding gene. Exhibits several functions, including double-stranded RNA binding activity; protein homodimerization activity; and transferrin receptor activity. Involved in several processes, including negative regulation of mitochondrial fusion; positive regulation of lymphocyte activation; and positive regulation of macromolecule metabolic process. Localizes to several cellular components, including HFE-transferrin receptor complex; clathrin-coated pit; and endosome. Implicated in breast carcinoma; immunodeficiency 46; multiple myeloma; and ovarian cancer. Biomarker of breast cancer; carcinoma (multiple); iron deficiency anemia; and prostate cancer.