



IVTScrip™ mRNA-Human ALDH1, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4799MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4799MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the ALDH1 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000689
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000689
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	ALDC; ALDH11; ALDH-E1
Description	This gene is a protein coding gene. Exhibits aldehyde dehydrogenase (NAD ⁺) activity. Involved in oxidation-reduction process. Localizes to cytosol. Biomarker of pancreatic cancer.