




IVTScrip™ mRNA-Human CD18, (Cap 1, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK4847MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK4847MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the CD18 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000211
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000211
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	LFA-1; MF17; LAD
Description	This gene is a protein coding gene. Exhibits several functions, including ICAM-3 receptor activity; heat shock protein binding activity; and protein kinase binding activity. Predicted to contribute to cargo receptor activity. Involved in several processes, including integrin-mediated signaling pathway; neutrophil chemotaxis; and positive regulation of neutrophil degranulation. Localizes to cell surface; integrin alphaL-beta2 complex; and plasma membrane raft. Colocalizes with cytoplasmic region. Implicated in coronary restenosis; leukocyte adhesion deficiency; and leukocyte adhesion deficiency 1. Biomarker of Vogt-Koyanagi-Harada disease and myocarditis.