



IVTScrip™ mRNA-Human MIXL1, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK5108MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5108MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the MIXL1 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001282402.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001282402.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MILD1; hMix; MGC138179; MIXL; MIX
Description	This gene is a protein coding gene. Exhibits RNA polymerase II repressing transcription factor binding activity and sequence-specific double-stranded DNA binding activity. Predicted to be involved in several processes, including endodermal cell differentiation; hematopoietic progenitor cell differentiation; and negative regulation of hematopoietic progenitor cell differentiation. Localizes to nucleus.