



IVTScrip™ mRNA-Human CD11a, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK5144MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5144MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CD11a protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_002209
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_002209
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CD11A; LFA1A
Description	This gene is a protein coding gene. Exhibits ICAM-3 receptor activity and cell adhesion molecule binding activity. Involved in several processes, including T cell activation via T cell receptor contact with antigen bound to MHC molecule on antigen presenting cell; heterophilic cell-cell adhesion via plasma membrane cell adhesion molecules; and memory T cell extravasation. Localizes to cell surface and integrin alphaL-beta2 complex.