



IVTScrip™ mRNA-Human CD220, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK5169MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5169MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CD220 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000208
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000208
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	IR; HHF5; CD220
Description	<p>This gene is a protein coding gene. Exhibits several functions, including insulin-like growth factor binding activity; peptide binding activity; and purine ribonucleoside triphosphate binding activity. Involved in several processes, including positive regulation of cellular metabolic process; positive regulation of glucose import; and protein phosphorylation. Localizes to caveola and insulin receptor complex. Colocalizes with nuclear envelope and nuclear lumen. Implicated in Donohue syndrome; diabetes mellitus (multiple); familial hyperinsulinemic hypoglycemia 5; and reproductive organ cancer (multiple). Biomarker of Alzheimer's disease; breast cancer; hepatocellular carcinoma; renal cell carcinoma; and seminoma.</p>