



## IVTScrip™ mRNA-Human CD99, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK5299MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK5299MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CD99 protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	5-Methyl-CTP & Pseudo-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001330312.19
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

### SPECIFICATIONS

<b>Modified bases</b>	5-Methyl-CTP & Pseudo-UTP
-----------------------	---------------------------

<b>5' Cap</b>	Cap 1
---------------	-------

<b>Species</b>	Human
----------------	-------

<b>RefSeq</b>	NM_001330312.19
---------------	-----------------

<b>Applications</b>	Gene therapy research
---------------------	-----------------------

<b>Format</b>	Powder
---------------	--------

<b>Quantity</b>	100 µg
-----------------	--------

<b>Purification</b>	oligo-dT affinity purification
---------------------	--------------------------------

## GENE INFORMATION

<b>Alternative Names</b>	MIC2Y; MSK5X
--------------------------	--------------

<b>Description</b>	This gene is a protein coding gene. Predicted to be involved in several processes, including T cell extravasation; myeloid dendritic cell activation; and positive regulation of neutrophil extravasation. Localizes to focal adhesion.
--------------------	---