




IVTScrip™ mRNA-Human CD69, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5462MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5462MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the CD69 protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001781
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001781
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	AIM; CLEC2C
Description	This gene is a protein coding gene. Exhibits identical protein binding activity. Predicted to be involved in cellular response to drug. Localizes to protein-containing complex.