



IVTScrip™ mRNA-Human IL17Ra, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5591MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5591MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the IL17Ra protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_014339
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_014339
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	IL17R; CD217
Description	This gene is a protein coding gene. Predicted to have interleukin-17 receptor activity. Involved in fibroblast activation; interleukin-17-mediated signaling pathway; and positive regulation of cytokine production. Predicted to localize to plasma membrane. Implicated in immunodeficiency 51.